

The product pluriPlix®

pluriPlix® is a rotation mixer which mixes all kinds of samples in a thorough and gentle way, and in a constant quality. This is essential for an optimal result of cell isolation. pluriPlix® does not have a drive with power connection of its own. Instead, it is driven by a standard magnetic stirrer and thus allowing for a variable speed of rotation. So, the pluriPlix® rotation mixer makes a flexible multi-purpose mixer of any standard magnetic stirrer.

pluriPlix® is available in a single-armed or a two-armed design. The magnetic lock facilitates changing the mixing options. The following options are available: over-head, tilting-rotation mixing and - when using the two-armed design - both options simultaneously.

pluriPlix® is suitable for all standard reaction tubes. The maximum mixing volume is 200 ml per arm (using 4x50 ml reaction tubes). It can be used within a temperature range of 0 to 45°C and also the rotation angle can be adjusted up to 10 degrees. No service required.